

*Trusted  
Your Partner in Healthcare*

# Kleen Mini

(Portable Suction Machine)



## FEATURES:

1. High negative pressure, large flux, low noise. It is a brand new mechanism design with suction bottle hidden inside and square negative pressure meter crust completely mad of plastic. This design not only reduces the working noise but also makes it more fashionable and precise.
2. It adopts advanced overflow-protection design, which can prevent pus and mucus from leaking into the machine.
3. Compact, lightweight, and easy to carry especially suitable for emergency and outside use.

## SPECIFICATIONS:

1. High vacuum, low flow
2. Power supply : AC220V $\pm$ 22V, 50Hz $\pm$ 1Hz
3. Input power : 90VA
4. Limit negative pressure : 0.075Mpa (760mmHg)
5. Negative pressure regulating range : 0.02MPa ~ limit negative pressure
6. Suction rate : 18L/min
7. Liquid holder : 1000ml/pc, 1pc
8. Sound level : 65dB(A)
9. Wight : 4.4kg
10. Dimension : 280x196x285(mm)
11. Reservoir Volume : 1000ml. 1 piece

- The apparatus is not suitable for use in the place with inflammable & explosive gas;
- Working system: short time running;
- Electric classification: Class II equipment, Type B applied part.

## NORMAL OPERATING CONDITIONS:

- Ambient temperature : +5°C +35°C
- Relative humidity : 30% ~ 80%
- Atmospheric pressure : 86 kPa ~106 kPa



## VIN BIOTECH SYSTEMS LIMITED

R.O.: No. 2, II Cross Street, M.E.S. Road, Tambaram (E), Chennai - 59.  
W: No. 2, U.V. Swaminathan Iyer Street, M.E.S. Road, Tambaram (E), Chennai - 600 059.  
Tel: +91 - 44 - 2239 0795, 6456 2736  
Email: info.vbsl@gmail.com, ceo.vinbiotech@gmail.com,  
Website: www.vinbiotech.com, Skype: shani.p7

