

*Trusted
Your Partner in Healthcare*

Kleen

(Suction Apparatus)



FEATURES:

Complete plastic panels, hand-switch and foot-switch for ease of operation.

Select negative pressure at will with no level for pressure adjustment.

Equipped with safety valve and overflow protection mechanism.

Safe and reliable, without pollution of oil.

Note: Not suitable for inflammable or explodable surroundings.

Max continuous working time is 30 minutes.

SPECIFICATIONS:

Max Vacuum	: 90kPa
Max airflow	: 20L/min
Adjustable vacuum range	: 20kPa ~ Max vacuum
Sound Level	: 65dB(A)
Suction bottle (glass)	: 2500ml/bottle, 2 Nos
Power supply	: AC230V \pm 10% 50Hz
Input power	: 90VA
Dimension	: 360 X 320X 480(mm ³)
Weight	: 15.5 kg
Electric classification	: I class B type

VIN BIOTECH SYSTEMS LIMITED

R.O.: No. 2, II Cross Street, M.E.S. Road, Tambaram (E), Chennai - 59.

W: No. 2, U.V. Swaminathan Iyer Street, M.E.S. Road, Tambaram (E), Chennai - 600 059.

Tel: +91 - 44 - 2239 0795, 6456 2736

Email: info.vbsl@gmail.com, ceo.vinbiotech@gmail.com,

Website: www.vinbiotech.com, Skype: shani.p7

